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Thank you for choosing the new Amplicord 55X electrolarynx, a simple and functional device to fit your lifestyle.

Contents of the Box [A]

- | | |
|--|------------------|
| ❶ Amplicord 55X | ❷ Security Strap |
| ❷ Power charger AC 220V 50Hz
DC 9V 30mA | ❸ Slip Case |
| ❸ 8.4V 250mAh rechargeable battery | ❹ Dial stopper |
| | ❺ User's Guide |

Product Description [B]

- | | |
|----------------------|------------------------|
| ❶ Device Head | ❷ Volume Dial V |
| ❷ Power button | ❸ Battery cover |
| ❸ Tone Dial T | ❹ Strap Hook |

Getting Started with Your Device [C]

1. Unscrew the battery cover.
2. Insert the batter, matching the + and - contacts.
3. If desired, attach the strap.
4. If desired, insert the dial stopper.
5. Screw the battery cover back on.

Recharging the Battery [D]

1. Unscrew the battery cover.
2. Remove the battery.
3. Insert the battery in the battery charger and plug it into a power outlet. It will take 8/10 hours to fully recharge the battery.
4. Remove the charged battery and insert it into the electrolarynx, aligning the contacts.
5. Screw the battery cover back on.

Instructions for Use [E]

1. Apply the device onto the softer part of the neck, with a light and even pressure.
2. Press the button to activate the electrolarynx and test it by vocalising the letters **A - E - I - O - U**.
3. *If you have inserted the dial stopper, slide it down to access the dials.*
4. Use the **T** dial to adjust to the desired tone.
5. Use the **V** dial to adjust to the desired volume.
6. *If you have inserted the dial stopper, slide it back up to lock the dials at the desired settings.*

Device Maintenance

Maintenance will be carried out by the dealer. The device may be cleaned as often as needed, using a piece of cloth and some rubbing alcohol (surgical spirit).

Warranty

The Amplicord 55X carries a 2-Year Limited Warranty (24 months) for each of its components, valid from the Date of Purchase. The device is guaranteed against manufacturing defects (defined as a defect in materials or workmanship). The warranty states that a free replacement and/or repair will be provided for any faulty components, during the applicable warranty period. This warranty will lose validity should the device show signs of manumission, or if the defect is a result of improper use or accidental drops or falls. The power charger and the rechargeable battery are excluded from this warranty.

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